


Multiple Tray options

The clinician can order the dreamTAP in a variety of trays.

Written By: Patrick Tessier



Step 1 — AMI Tray systems



- **AccuTherm trays:** The new option. Kind of a mashup of TL and ThermAcryl.
 - 1.0 mm of a PET-G shell, bonded with a thin layer of ThermAcryl.
 - Can be remolded for tighter or looser fit. Can even add ThermAcryl to increase retention.
- **Triple Laminate (TL) trays:** These are the most popular trays. They feature a soft urethane liner and a poly-carbonate shell.
 - Great for fast delivery. They literally snap on
- **ThermAcryl trays:** These are popular with seasoned DSM clinicians. They feature a heat moldable liner called ThermAcryl.
 - Excellent choice for short teeth with retention issues.
 - Also great if you are going to modify for future restorations.