

ThermAcryl - White

Written By: Patrick Tessier

Comments: 4
Favorites: 0
Completions: 1



🔒 Difficulty 🟡 🟡 ⚪ Moderate
☰ Steps 1
🕒 Time Required Suggest a time ??
🔖 Sections 1 ▾
⚠️ Flags 0 ▾

Introduction

- This is our regular ThermAcryl
- Turns clear and soft at about 160 degrees F
 - Little bit sticky
 - 2% shrink factor

Step 1 ThermAcryl - White



- ThermAcryl comes in 2 packages
 - 2 Oz bags
 - 1 Lbs 10oz tins
- Uses
 - Add on to myTAP trays
 - Add on to ThermAcryl lined trays
 - Add on to AccuTherm trays

[Add a comment](#)

FINISH LINE

[✔ Complete Guide](#)

One other person completed this guide.

<p>Author</p> <div style="display: flex; align-items: center;"> <div> <p>Patrick Tessier</p> <p>Member since: 08/09/2017</p> <p>56 Guides authored</p> </div> </div>	<p>Team</p> <div style="display: flex; align-items: center;"> <div> <p>AMI</p> <p>Community</p> <p>11 Members</p> <p>56 Guides authored</p> </div> </div>
--	---

4 COMMENTS [Add a comment](#)

- 💬 i would like to order thermacryl - white, please advise on how to proceed. i do not see any contact information to call and ask for a quote.

[Ivan Vargas](#) - January 28
- 💬 Howdy Ivan,
You can call our main number to order: 866-264-7667
Fees: \$46.50 for a can (26 oz) or \$4.75 for a 2 oz bag.
Patrick
patrick@amisleep.com

[Patrick Tessier](#) - January 28
- 💬 Please give me a call at 520-533-5656 thank you.

[Ivan Vargas](#) - January 28
- 💬 is it possible if you can email me a quote to my email? Please include shipping cost. I would like 4 bags of Thermacryl white. Below is my information:
ivan.w.vargas.civ@mail.mil
Phone: 520-533-5656
Thank you!

[Ivan Vargas](#) - January 28

Add Comment

B I 🔗 ” | ↶ ↷

[Post comment](#)