

Adjusting a dreamTAP ThermAcryl

Written By: Patrick Tessier



INTRODUCTION

Since a ThermAcryl lined TAP appliance has the final fit in the dental practice, any adjustments are done following the same procedures.



PARTS:

- [ThermAcryl Material](#) (1)
-

Step 1 — Adjusting for poor fit



- If the patient is not satisfied with the fit of the appliance, the adjustment procedure is the same as the fit procedure. [Delivering a dreamTAP ThermAcryl guide.](#)

Step 2 — Complete relining of device



- If the adjustments are major, as in the case of a new restoration or major fit issues, you can reline the appliance and basically start from the beginning. Here is the [Relining with ThermAcryl Guide](#).