

NEW AM Aligner™



KEY FEATURES

- Handle (tab) allows for easier initial fitting and daily use
- Duplicates patient's bite before starting therapy
- Matches habitual occlusion
- Shorter heating time and fitting time compared to the classic AM Aligner



tapintosleep.com




contactami@amisleep.com



1-866-264-7667

WHAT IS THE AM ALIGNER?

Mandibular protrusion can move the jaw forward permanently over time if the patient doesn't make sure that normal occlusion is reestablished each day.




Included with every TAP appliance, the **AM Aligner** is a unique thermoplastic material that takes a registration of the patient's normal bite before treatment. It is to be used for 5-10 minutes each morning after use of an oral appliance to exercise the jaw back to its usual position.

HOW TO FIT

Fitting the new AM Aligner takes less than 10 minutes.

Watch the below video or refer to the AM Aligner instructions for use for fitting instructions.



Scan to watch fitting video