

ORDERING A CUSTOM TAP™ APPLIANCE

PRINT A LAB SLIP & SHIP, OR SUBMIT DIGITAL SCANS

When you are ready to order your first case, print out an [Airway Labs Lab Slip](#) and ship your models using our [pre-paid UPS shipping labels](#), or [submit your digital scans](#) using our online form. Go to www.airwaylabs.com under “How To Order” to find the most up-to-date digital lab slip and a list of frequently asked questions.

When you complete an Airway Labs Prescription form and send it with the patient’s models or Poly Vinyl impressions, we will call to inform you that we have received your case and obtain your contact information to **create an account for you**. Your account number will appear on your invoice. When we return the appliance to you, we will include our lab slips and a UPS label (for future cases or returns).

CHOOSE HARDWARE AND LINER

flexTAP

- Vertex Technology “opens the box” with less protrusion
- More tongue space due to less hardware in the mouth
- Standard Mouth Shield nasal breather trainer slips onto post
- 14 mm range of self-titratable adjustment
- No metal parts
- New ThermAcryl liner is thinner and easier to work with



dreamTAP

- Nickel-free, allergen-free metal
- 15 mm range of advancement (3 hook sizes)
- Custom-formed soft liners for comfort and superior retention (AccuTherm or Triple Laminate)



TAP 3

- Single hook allows for 7mm range of advancement
- Stainless steel parts
- Custom-formed soft liners for comfort and superior retention (AccuTherm or Triple Laminate)

