




**Morning exercise tool to return the bite back to its usual position**



### KEY FEATURES

The AM Aligner® is a unique thermoplastic material that takes a registration of the patient's natural bite before treatment.

- Tab allows for easier initial fitting and daily use
- ThermAcryl® material is thinner, comfortable, resilient, and flexible
- Shorter heating time and fitting time compared to the classic AM Aligner
- Anterior teeth guide helps patients bite down properly

### HOW TO FIT

Fitting the AM Aligner takes less than 10 minutes.

New white color turns clear when ready to fit.

Scan the QR code to watch the fitting video.



### WHY USE THE AM ALIGNER?

All mandibular advancement devices can move the jaw forward permanently over time if the patient doesn't ensure that natural occlusion is reestablished each day.

The AM Aligner is to be used for 5-10 minutes each morning after use of an oral appliance to exercise the mandible back to its usual position.



Standard with every TAP® appliance

### WANT TO LEARN MORE?

Email [contactami@amisleep.com](mailto:contactami@amisleep.com) for more information.

Visit our website [www.tapintosleep.com](http://www.tapintosleep.com).

